

# Analyse case7 dicompyler-core 28-11-2016

November 28, 2016

Metrics show high PTV dose and low brain stem dose in this case

```
In [1]: %matplotlib inline
import os
import numpy as np
import dicom as dicom
from dicompylercore import dicomparser, dvh, dvhcalc
import matplotlib.pyplot as plt
import urllib.request
import os.path
from ipywidgets import FloatSlider, interactive, interact
from scipy.misc import imread
from scipy import interpolate
from pylinac.core.image import ArrayImage
import copy
import scipy.ndimage.filters as spf
```

## 1 Helpers

```
In [3]: def f_info(file_):      # pass a dicomparser object, e.g. AAA
    print('****')
    print('position = ' + str( file_.GetDoseData()['position'] ))
    print('dosegridscaling = ' + str( file_.GetDoseData()['dosegridscaling'] ))
    print('frames = ' + str( file_.GetDoseData()['frames'] ))
    print('rows = ' + str( file_.GetDoseData()['rows'] ))
    print('columns = ' + str( file_.GetDoseData()['columns'] ))
    print('pixelspacing = ' + str( file_.GetDoseData()['pixelspacing'] ))
    #print('GridFrameOffsetVector = ' + str( file_.ds.GridFrameOffsetVector ))
    print('dosemax (pix) = ' + str( file_.GetDoseData()['dosemax'] ))
    print('dosemax (Gy) = ' + str( file_.GetDoseData()['dosemax'] * file_.C ))
    #print('samplesperpixel = ' + str( file_.GetDoseData()['samplesperpixel'] ))

In [4]: def get_dose_array(rtdose):      # e.g. rtdose = AAA dicomparser object, return
    x = rtdose.GetDoseData()['rows']
    y = rtdose.GetDoseData()['columns']
    z = rtdose.GetDoseData()['frames']
    rtdose_array = np.zeros((x, y, z))    # create 3d array for data
```

```

planes = (np.array(rtdose.ds.GridFrameOffsetVector) * rtdose.ds.ImageOn
for i, plane in enumerate(planes): # put dose data in 3d array
    rtdose_array[:, :, i] = rtdose.GetDoseGrid(plane)
rtdose_array = rtdose_array * rtdose.GetDoseData()['dosegridscaling']
return(rtdose_array)

# plt.imshow(get_dose_array(AAA)[:,:,: 50]); # exampleusage

```

## 1.1 Analysis

```

In [5]: Structures_set = dicomparser.DicomParser('Case7_structures.dcm')
        AAA = dicomparser.DicomParser('Case7_dose_AAA.dcm') # a dicomparser object
        Dm = dicomparser.DicomParser('Case7_dose_Dm.dcm')

```

```

In [6]: f_info(AAA)

```

```

****
position = ['-244.4338517', '-178.5440886', '-216.25']
dosegridscaling = 6.501e-5
frames = 135
rows = 107
columns = 192
pixelspacing = ['2.5', '2.5']
dosemax (pix) = 1066131.0
dosemax (Gy) = 69.30917631

```

```

In [7]: f_info(Dm)

```

```

****
position = ['-245.6838517', '-179.7940886', '-216.25']
dosegridscaling = 7.8024484e-5
frames = 135
rows = 108
columns = 193
pixelspacing = ['2.5', '2.5']
dosemax (pix) = 942825.0
dosemax (Gy) = 73.5634341273

```

### 1.1.1 Get dose grids, crop array to match

```

In [8]: AAA_dose = get_dose_array(AAA)
        print(AAA_dose.shape)
        x = AAA_dose.shape[0]
        y = AAA_dose.shape[1]
        print('***')

        Dm_dose = get_dose_array(Dm)

```

```

print(Dm_dose.shape)
Dm_dose = Dm_dose[:,y, :] # put on same grid
print(Dm_dose.shape)

(107, 192, 135)
***
(108, 193, 135)
(107, 192, 135)

In [9]: y = FloatSlider(
        value=75,
        min=0,
        max=135,
        step=1,
        description='Slice Position in numpy:',
    )

def showdoseboth(z):

    AAA_display = AAA_dose[:, :, z]
    Dm_display = Dm_dose[:, :, z]
    cmin = min(AAA_display.min(), Dm_display.min()) # use identical colour
    cmax = max(AAA_display.max(), Dm_display.max()) # use identical colour

    f, (ax1, ax2, ax3) = plt.subplots(1, 3, figsize=(12, 18))

    ax1.imshow(AAA_display, vmin = cmin, vmax = cmax) # , vmin=0, vmax=cmax
    ax1.set_title('AAA, mean = ' + str(np.mean(AAA_display)))


    ax2.imshow(Dm_display, vmin = cmin, vmax = cmax) # , vmin=0, vmax=cmax
    ax2.set_title('Dm, mean = ' + str(np.mean(Dm_display)))

    # Difference analysis
    difference_display = (Dm_display - AAA_display)
    ax3.imshow(difference_display, cmap='hot', vmin = 0.0, vmax = 10.0) #
    ax3.set_title('(Dm - AAA), mean = ' + str(np.mean(difference_display)))

    plt.show()

interactive(showdoseboth, z=75)

```



## 1.2 DVH analysis

```
In [10]: structures = Structures_set.GetStructures() # returns a dict of structures
```

```
for i, structure in enumerate(structures):
    print(str(i+1) + ' ' + structures[structure]['name'])
```

```
1 GTV_P8
2 GTVn_P2
3 GTV + 15mm_P1
4 GTVn + 10mm_P1
5 CTV65_P12
6 CTV54_P8
7 L Parotid_P14
8 R Parotid_P14
9 L Lens_P12
10 R Lens_P12
11 L Orbit_P12
12 R Orbit_P12
13 Brain Stem_P14
14 Spinal Cord_P14
15 BS + 3mm_P9
16 BS + 5mm_P13
17 Body_P14
18 Lt Par edit_P14
19 PTV1_P13
20 PTV1 edit_P13
21 PTV2_P13
22 PTV2 edit_P13
23 Pseudo Midline14
24 Pseudo Oral Ca14
25 Pseudo PostBrai2
26 Pseudo Post_P12
27 Rt Par edit_P14
28 SC + 3mm_P14
29 SC + 5mm_P14
30 Pseudo PTV1_P14
31 Pseudo PTV2_P14
```

```

32 Left Couch Bar_6
33 Right Couch Bar6
34 Artifact_P2
35 Annulus1_P3
36 Dose 6175[cGy]_4
37 Dose 5130[cGy]_4
38 missing1_P1
39 missing2_P1

```

```

In [12]: def compare_dvh(key=1):
          structure = Structures_set.GetStructures()[key]
          AAA_ = dvh.DVH.from_dicom_dvh(AAA.ds, key, name=structure['name'] + 'AAA')
          Dm_ = dvh.DVH.from_dicom_dvh(Dm.ds, key, name=structure['name'] + 'Dm')
          AAA_.compare(Dm_)


```

```

In [13]: compare_dvh(key=30)


```

Structure:	Pseudo PTV1_P14AAA	Pseudo PTV1_P14Dm	Rel Diff	Abs diff
-----				
DVH Type:	cumulative, abs dose: Gy, abs volume: cm3			
Volume:	504.47 cm3	504.47 cm3	+0.00 %	+0.00
Max:	69.31 Gy	73.57 Gy	+6.15 %	+4.26
Min:	59.58 Gy	57.79 Gy	-3.00 %	-1.79
Mean:	64.99 Gy	67.37 Gy	+3.67 %	+2.38
D100:	0.00 Gy	0.00 Gy	+0.00 %	+0.00
D98:	62.59 Gy	64.43 Gy	+2.94 %	+1.84
D95:	63.18 Gy	65.06 Gy	+2.98 %	+1.88
D2cc:	67.68 Gy	71.90 Gy	+6.24 %	+4.22



```
In [14]: compare_dvh(key=13)
```

Structure:	Brain Stem_P14AAA	Brain Stem_P14Dm	Rel Diff	Abs diff
-----				
DVH Type:	cumulative,	abs dose: Gy,	abs volume: cm3	
Volume:	27.62 cm3	27.62 cm3	+0.00 %	+0.00
Max:	48.97 Gy	50.64 Gy	+3.41 %	+1.67
Min:	3.45 Gy	2.66 Gy	-22.90 %	-0.79
Mean:	15.97 Gy	15.37 Gy	-3.76 %	-0.60
D100:	0.00 Gy	0.00 Gy	+0.00 %	+0.00
D98:	3.71 Gy	2.82 Gy	-23.99 %	-0.89
D95:	3.90 Gy	2.96 Gy	-24.10 %	-0.94
D2cc:	41.52 Gy	42.79 Gy	+3.06 %	+1.27



In [ ]: